



## CERTIFICATE OF ANALYSIS

### Apolipoprotein A-II

#### Product Specifications

---

<b>Product code:</b>	<b>P230-0</b>
<b>Batch No:</b>	65-138
<b>Product presentation:</b>	Supplied as a 0.2 $\mu$ m filtered frozen solution in 48mM sodium phosphate buffer pH 7.2 containing 0.145M NaCl, 3mM Benzamidine Hydrochloride Monohydrate and 0.2% sodium azide.
<b>Purity:</b>	Nominally greater than 99% of the protein is Apolipoprotein A-II.
<b>Starting material:</b>	Human plasma. Starting material donor tested and found negative for HIV I & II antibodies, Hepatitis B surface antigen, and Hepatitis C antibodies.
<b>Storage:</b>	Store below $-15^{\circ}\text{C}$ .
<b>Expiry Date:</b>	22/April/2011

#### Product Analysis

---

<b>SDS-Polyacrylamide Gel Electrophoresis:</b>	To show one band corresponding to Apolipoprotein A-II molecular weight.
<b>Total protein concentration as determined by Behring Nephelometer-System:</b>	10 mg/ml